


Tracheostomy



ID	OD	Un	PR	PL	Adj.
3.0	5.0	0			
3.5	5.6	0			0
4.0	6.1	0	0		0
4.5	6.6	0	0		0
5.0	7.3	0	0	0	0
5.5	8.0	0	0	0	
6.0	8.7	0	0	0	0
6.5	9.3	0	0	0	
7.0	10.0	0	0	0	0
7.5	10.6	0	0	0	
8.0	11.2	0	0	0	0
8.5	11.7		0	0	
9.0	12.3		0		

*Tracheostomy Type

Un	UnCuff
PR	Preformed Cuff
PL	Plain Cuff
Adj	Adjustable