


# SILICONE STOMACH (GASTRIC) TUBE

Fortune® Silicone Stomach (Gastric) Tube is designed for nutrition supplement to stomach and may be recommended for various purposes: for patients who cannot take food or swallow, take enough food by mouth to keep nutrition, congenital defects of mouth, esophagus, or stomach.

## Fortune® Silicone Stomach (Gastric) Tube Features:

- Medical grade silicone for superior biocompatibility.
- Proper softness and excellent flexibility for enhanced patient comfort.
- Tapered tip for easy insertion.
- Depth marking for position confirmation.



## SPECIFICATIONS

### Stomach (Gastric) Tube:

REF. No.	Size	Length	Description
2020-0004	4FR	400 mm	
2020-0005	5FR		- Attached plug
2020-0006	6 FR		- Depth marking - X-ray opaque line
2020-0008	8 FR	800 mm	
2020-0010	10 FR		

REF. No.	Size	Length	Description
2020-0012	12 FR	1200 mm	
2020-0014	14 FR		- Attached plug
2020-0016	16 FR		- Depth marking - X-ray opaque line
2020-0018	18 FR		
2020-0020	20 FR		



2020-00XX